


# NECK JUNCTIONAL HEMORRHAGE CONTROL



**CONSIDER** body substance isolation.  
**NOTE:** If a Combat Lifesaver is available, direct them to assist.

**01 EXPOSE** the injury and assess the bleeding source.  
**NOTE:** Clothing may need to be cut away to properly expose the injury.



**02 APPLY DIRECT PRESSURE** to bleeding source.

**03 REMOVE** hemostatic dressing from casualty's JFAK.

**NOTE:** If hemostatic dressing is not available, use gauze or clean, dry cloth material.



**04 PACK** wound tightly with hemostatic dressing until the wound cavity is filled. Continue to apply pressure. Finish packing within 90 sec.

**NOTE:** More than one dressing/gauze may be required to stop blood flow.

**CAUTION:** If a penetrating object is lodged in the casualty's body, bandage it in place. **DO NOT** remove the object.




**05 ENSURE** hemostatic dressing or gauze extends 1-2" above the skin.

**NOTE:** If the hemostatic dressing or gauze does not extend 1-2" above the skin, place additional hemostatic dressing or gauze.



**06** After packing, continue to apply firm, manual pressure until the bleeding stops. **HOLD** pressure for a minimum of 3 min.



**07 REASSESS** to ensure bleeding has been controlled while maintaining pressure.


**NOTE:** Ensure blood is not seeping through or around the hemostatic dressing. If bleeding has stopped, leave the hemostatic dressing in place while still maintaining pressure.

**IF BLEEDING HAS NOT BEEN CONTROLLED:**

**08a** If packed with hemostatic dressing, remove prior packing material and repack starting at **STEP 3**.

**OR**

**08b** If packed with gauze, apply additional gauze/materials and pressure (for 3 min) until bleeding stops.



**09 PLACE** a 6" elastic bandage over the hemostatic dressing, leaving enough tail to tie it into a knot.



**10** While maintaining pressure, **WRAP** bandage (no less than 1½ times) over the packing material covering it completely.



**WRAP** diagonally across the chest under the opposite arm (axilla) around the back to the neck and back over the wound.

*Continued on the next page...*


# NECK JUNCTIONAL HEMORRHAGE CONTROL

*Continued...*




**11 WRAP** around neck and under the arm (on alternating sides of the tail, while maintaining tension/pressure), pulling elastic bandage tightly for pressure, covering the packing material.

**NOTE:** Do not use pressure bar on the neck.



**12 SECURE** the dressing by tying a non-slip knot with end of elastic bandage and its tail.




**13 SECURE** elastic bandage tails with 3" tape, wrapping a minimum of 1½ times around the knot.



**14 SWATHE** the upper arm (of the uninjured side) to the chest using a bandage.



**15 CONTINUE TO ASSESS** wound for further bleeding.



**16 DOCUMENT** all findings and treatments on a DD Form 1380 TCCC Casualty Card and attach it to the casualty.